

STEAMJET

Klasse B
AUTOCLAVE

15
min

Total process duration for 134°C process setting is just 15 minutes. Ultra short sterilization time allows preparing tools when they are needed.




2 YEARS
WARRANTY



The large touchscreen facilitates intuitive use. In addition to basic functions, it is possible to set up your own sterilization cycle.

The process is presented in a graphic form giving information on remaining time, temperature and pressure.



STEAMJET



2 YEARS
WARRANTY


All STEAMJET autoclaves are covered by a full two year warranty.

Program	134°C	121°C
Tools	wrapped and unwrapped	
Sterilization time	4 min	15 min
Drying time	3 min	5 min
Total process duration	15 min	30 min

Klasse B
AUTOCLAVE

TECHNICAL DATA

Power supply	230V/50Hz
Power output	3500 W
Maximum current consumption	15 A
Operating pressure	2,1 Bar
Maximum pressure	2,3 Bar
Maximum process temperature	137°C
Volume of the process chamber	2,7 l
Process chamber dimensions (LxWxH)	292 x 192 x 45mm
External dimensions (LxWxH)	561 x 252 x 162mm
Weight	15 kg
IP Code	IP20
Process data logging	USB Drive



Enbio Group AG
Zugerstrasse 76 B, 6340 Baar, Switzerland
T. +49 30 863 20 192, info@enbio-group.com
www.enbio-group.com

ver: 06.12.2018

enbio

STEAMJET

Klasse B
AUTOCLAVE

Autoclave for tools

100% protection against infections



the
highest
quality

www.enbio-group.com

enbio

STEAMJET

Compact

Ultraquiet

Rapid

Intuitive

Sterilisation on demand

Compact design, small size, extremely quiet work and intuitive operation of two programs. Autoclave STEAMJET is the best solution for beauty salons. Instant sterilization process makes each instrument available just in time.

In every beauty salon there is a risk of dangerous viral infections, e.g. jaundice, scabies and even HIV. Thanks to the Steamjet, tools are always sterile and customers can feel safely. Dangerous viruses have no chance.

Lowest price: 985 Euro net

Made in Switzerland

The highest quality of sterilisation

