

Enbio S - Autoclaves for sterilizing tools

Sterile and safe practice



**TOP
QUALITY**





2 YEARS WARRANTY

100% PROTECTION

MEDICINE PRODUCT

UNBEATABLE PRICE

Program

	134°C FAST*	134°C	121°C
Sterilization time [min]	3,5	4	15
Drying time [min]	-	3	5
Total process duration [min]	7	15	30

*The FAST 134 (Class S Sterilization) program is for flat, without porous tools which are sterilized only unpackaged and without any textile.



EASY TO USE



SUPERFAST



COMPACT

Technical data

Power supply	230V/50Hz
Power output	3250 W
Maximum power consumption	15 A
Operating pressure	2,1 Bar
Maximum pressure	2,3 Bar
Maximum process temperature	137°C
Chamber volume	2,7 l
Max. number of packs (about)	8
Load size	0,5 kg
Weight of the device	15 kg
IP code	IP20
Archiving the process data	USB Drive
Connection of the water canisters	quick connector
Demineralized or distilled water	external



Sales in Europe

Enbio Group GmbH
Basler Str. 115
79115 Freiburg, Germany



Headquarters in Switzerland

Enbio Group AG
Eichengasse 3
4702 Oensingen, Switzerland

T. +49 761 88793141
info@enbio-group.com
www.enbio-group.com

ver: 17.07.2019

